

HYBRID FIXATOR
LARGE BAR FIXATOR
RAIL FIXATOR

RING FIXATOR

Standard, Small and Ankle with Hybrid and Ring Fixator compatibility



Combining Strength, Versatility and Flexibility the UniX Unilateral and Ring Fixator System

- Standard Unilateral Fixator
- Ankle and T-Clamp options
- Mechanical locking of articulating components allows for non-slip locking
- Modular system - maximum flexibility

Part of the Ring Fixator Family, allowing easy conversion to a hybrid or complete ring fixation system

*MR
compatible*

The UniX Unilateral Fixator

- Serrated locking connectors provide a non-slip locking mechanism.
- Connectors each provide 120° of articulation.
- Central component provides 360° rotation.
- Up to 0-5 cm of compression and/or distraction in each clamp, for a total of 10cm of compression/distraction capability.
- Modular system maximum flexibility.



UniX Ankle Clamp

- Provides acute and/or long-term treatment option.
- Compatible with percutaneous reduction and fixation techniques.
- Proven C-clamp design allows imaging compatibility.
- No templates required.
- Serrated connectors provide a non-slip locking mechanism.
- Modular system.
- Maximum flexibility.



Telescopic Strut

- Provide to convert the frame to the Ilizarov type ring fixator easily
- Gradually adjustable telescopic design.

UniX Hybrid Clamp

- Full, 2/3 and Foot Rings.
- Aluminium and Carbon Fiber constructs.
- Fully convertible to the Ring Fixator system.



Advanced Screw Clamp

- Provide desired screw position on the ring
- Fully adjustable screw height angles.



Bolt, Nut, Washer

- Titanium bolt, nut and washer provide lightweight frame and improve MRI compatibility.



ORDERING INFORMATION

Catalog No	Description
Ring Fixation System	
50-1012-10A	Dual Hole Full Ring, 100mm inner diameter, Aluminium
50-1012-12A	Dual Hole Full Ring, 120mm inner diameter, Aluminium
50-1012-13A	Dual Hole Full Ring, 135mm inner diameter, Aluminium
50-1012-15A	Dual Hole Full Ring, 150mm inner diameter, Aluminium
50-1012-16A	Dual Hole Full Ring, 165mm inner diameter, Aluminium
50-1012-18A	Dual Hole Full Ring, 180mm inner diameter, Aluminium
50-1012-19A	Dual Hole Full Ring, 195mm inner diameter, Aluminium
50-1012-21A	Dual Hole Full Ring, 210mm inner diameter, Aluminium
50-1012-22A	Dual Hole Full Ring, 225mm inner diameter, Aluminium
50-1017-12A	Dual Hole 2/3 Ring, 120mm inner diameter, Aluminium
50-1017-13A	Dual Hole 2/3 Ring, 135mm inner diameter, Aluminium
50-1017-15A	Dual Hole 2/3 Ring, 150mm inner diameter, Aluminium
50-1017-16A	Dual Hole 2/3 Ring, 165mm inner diameter, Aluminium
50-1017-18A	Dual Hole 2/3 Ring, 180mm inner diameter, Aluminium
50-1017-19A	Dual Hole 2/3 Ring, 195mm inner diameter, Aluminium
50-1017-21A	Dual Hole 2/3 Ring, 210mm inner diameter, Aluminium
50-1017-22A	Dual Hole 2/3 Ring, 225mm inner diameter, Aluminium
50-1012-12C	Dual Hole Full Ring, 120mm inner diameter, Carbon Composite
50-1012-13C	Dual Hole Full Ring, 135mm inner diameter, Carbon Composite
50-1012-15C	Dual Hole Full Ring, 150mm inner diameter, Carbon Composite
50-1012-16C	Dual Hole Full Ring, 165mm inner diameter, Carbon Composite
50-1012-18C	Dual Hole Full Ring, 180mm inner diameter, Carbon Composite
50-1012-19C	Dual Hole Full Ring, 195mm inner diameter, Carbon Composite
50-1012-21C	Dual Hole Full Ring, 210mm inner diameter, Carbon Composite
50-1012-22C	Dual Hole Full Ring, 225mm inner diameter, Carbon Composite
50-1031-00	Dual Joint Strut, Extra Extra Short (XXS)
50-1031-01	Dual Joint Strut, Extra Short (XS)
50-1031-02	Dual Joint Strut, Short (S)
50-1031-03	Dual Joint Strut, Medium (M)
50-1031-04	Dual Joint Strut, Long (L)
50-1031-05	Dual Joint Strut, Extra Long (XL)
50-1032-00	Dual Joint Express Strut, Extra Extra Short (XXS)
50-1032-01	Dual Joint Express Strut, Extra Short (XS)
50-1032-02	Dual Joint Express Strut, Short (S)
50-1032-03	Dual Joint Express Strut, Medium (M)
50-1032-04	Dual Joint Express Strut, Long (L)
50-1041-30	Advanced Pin Clamp, Short
50-1041-00	Pin Clamp
50-1042-30	Advanced Screw Clamp, Short
50-1042-40	Advanced Screw Clamp, Long
50-1042-00	Screw Clamp
50-1051-00	Washer
50-1052-01	Standard Nut
50-1052-02	NyLock Nut
50-1053-10	Connection Bolt, XShort
50-1053-12	Connection Bolt, Short
50-1053-16	Connection Bolt, Medium
50-1053-20	Connection Bolt, Long
50-1053-30	Connection Bolt, X-Long
50-1054-20	Dual Rod Connector, Short
50-1054-40	Dual Rod Connector, Medium
50-1054-60	Dual Rod Connector, Long
50-1055-45	Twisted Plates, Short
50-1055-65	Twisted Plates, Medium
50-1055-85	Twisted Plates, Long
50-1058-06	Set Screw, 6mm
50-1058-04	Set Screw, 4mm
50-1041-90	Pin Fixation Bolt
50-1011-13	Half Ring, 135mm
50-1011-15	Half Ring, 150mm
50-1011-16	Half Ring, 165mm
50-1011-18	Half Ring, 180mm
50-1011-19	Half Ring, 195mm
50-1011-21	Half Ring, 210mm
50-1011-22	Half Ring, 225mm
50-1013-14	Femoral Arch, 90°, Small
50-1013-18	Femoral Arch, 90°, Large
50-1014-14	Femoral Arch, 120°, Small
50-1014-18	Femoral Arch, 120°, Large
50-1056-01	Cube, 1 hole
50-1056-02	Cube, 2 hole
50-1056-03	Cube, 3 hole
50-1056-04	Cube, 4 hole
50-1056-05	Cube, 5 hole
50-1056-60	Screw Sleeve for Cube, 6mm
50-1015-01	Foot Ring, Small, Aluminium
50-1015-02	Foot Ring, Large, Aluminium
50-1057-60	Threaded Rods, 60mm length
50-1057-80	Threaded Rods, 80mm length
50-1057-100	Threaded Rods, 100mm length

Catalog No	Description
50-1057-120	Threaded Rods, 120mm length
50-1057-150	Threaded Rods, 150mm length
50-1057-200	Threaded Rods, 200mm length
50-1057-300	Threaded Rods, 300mm length
UniX Unilateral Fixator	
50-5011-00	Unix Multijoint External Fixator, XS Model
50-5011-01	Unix Multijoint External Fixator C/D Mechanism, XS Model
50-5012-00	Unix Multijoint External Fixator, Small Model
50-5012-01	Unix Multijoint External Fixator C/D Mechanism, Small Model
50-5013-00	Unix Multijoint External Fixator, Large Model
50-5013-01	Unix Multijoint External Fixator C/D Mechanism, Large Model
50-5021-00	Unix Large Convergent T-Clamp, Offsett Screws
50-5022-00	Unix Large Convergent T-Clamp, Parallel Screws
50-5023-00	Unix Large T-Clamp, Standard
50-5021-01	Unix Small Convergent T-Clamp, Offsett Screws
50-5022-01	Unix Small Convergent T-Clamp, Parallel Screws
50-5023-01	Unix Small T-Clamp, Standard
50-5024-00	Unix Ankle C Clamp
50-5025-00	Unix Hip Distraction Kit
50-5026-02	Unix Long Clamps, 8cm, Right
50-5026-01	Unix Long Clamps, 8cm, Left
50-5027-00	Unix Anterior Hinge Clamp
50-5028-00	Unix Ring Adapter
50-5029-00	Unix Supplemental Rod Clamp
50-5028-15	Unix Supplemental Rod, Short
50-5028-30	Unix Supplemental Rod, Long
External Fixation Wires and Screws	
50-1061-18	Wire - 1.8mm
50-1061-20	Wire - 2.0mm
50-1061-18S	Wire - 1.8mm
50-1061-20S	Wire - 2.0mm
50-1063-18	Wire - 1.8mm, Trocar tip
50-1064-20	Wire - 2.0mm, Trocar tip
50-1063-18S	Wire - 1.8mm, Trocar tip
50-1064-20S	Wire - 2.0mm, Trocar tip
50-1062-18	Wire with Stopper- 1.8mm
50-1062-20	Wire with Stopper- 2.0mm
50-1062-18S	Wire with Stopper- 1.8mm
50-1062-20S	Wire with Stopper- 2.0mm
50-1065-220	Threaded Tip Wire, 2mm diameter, 20mm thread, 300mm length
50-1065-320	Threaded Tip Wire, 2mm diameter, 30mm thread, 300mm length
50-1065-420	Threaded Tip Wire, 2mm diameter, 40mm thread, 300mm length
50-1065-520	Threaded Tip Wire, 2mm diameter, 50mm thread, 300mm length
50-1065-225	Threaded Tip Wire, 2.5mm diameter, 20mm thread, 300mm length
50-1065-325	Threaded Tip Wire, 2.5mm diameter, 30mm thread, 300mm length
50-1065-425	Threaded Tip Wire, 2.5mm diameter, 40mm thread, 300mm length
50-1065-525	Threaded Tip Wire, 2.5mm diameter, 50mm thread, 300mm length
50-1074-122	Bone Screw, 4.5mm/120mm, 20mm thread, Trocar Point
50-1074-123	Bone Screw, 4.5mm/120mm, 30mm thread, Trocar Point
50-1074-124	Bone Screw, 4.5mm/120mm, 40mm thread, Trocar Point
50-1074-153	Bone Screw, 4.5mm/150mm, 30mm thread, Trocar Point
50-1074-154	Bone Screw, 4.5mm/150mm, 40mm thread, Trocar Point
50-1076-103	Bone Screw, 6.0mm/100mm, 30mm thread, Trocar Point
50-1076-113	Bone Screw, 6.0mm/110mm, 30mm thread, Trocar Point
50-1076-153	Bone Screw, 6.0mm/150mm, 30mm thread, Trocar Point
50-1076-154	Bone Screw, 6.0mm/150mm, 40mm thread, Trocar Point
50-1076-183	Bone Screw, 6.0mm/180mm, 30mm thread, Trocar Point
50-1076-184	Bone Screw, 6.0mm/180mm, 40mm thread, Trocar Point
50-1076-203	Bone Screw, 6.0mm/200mm, 30mm thread, Trocar Point
50-1076-204	Bone Screw, 6.0mm/200mm, 40mm thread, Trocar Point
50-1076-205	Bone Screw, 6.0mm/200mm, 50mm thread, Trocar Point
50-1076-223	Bone Screw, 6.0mm/220mm, 30mm thread, Trocar Point
50-1076-224	Bone Screw, 6.0mm/220mm, 40mm thread, Trocar Point
50-1076-225	Bone Screw, 6.0mm/220mm, 50mm thread, Trocar Point
50-1076-253	Bone Screw, 6.0mm/250mm, 30mm thread, Trocar Point
50-1076-254	Bone Screw, 6.0mm/250mm, 40mm thread, Trocar Point
50-1074-102B	Bone Screw, 4.5mm/100mm, 20mm thread, Blunt Point
50-1074-104B	Bone Screw, 4.5mm/100mm, 40mm thread, Blunt Point
50-1074-122B	Bone Screw, 4.5mm/120mm, 20mm thread, Blunt Point
50-1074-123B	Bone Screw, 4.5mm/120mm, 30mm thread, Blunt Point
50-1074-124B	Bone Screw, 4.5mm/120mm, 40mm thread, Blunt Point
50-1074-153B	Bone Screw, 4.5mm/150mm, 30mm thread, Blunt Point
50-1074-154B	Bone Screw, 4.5mm/150mm, 40mm thread, Blunt Point
50-1076-103B	Bone Screw, 6.0mm/100mm, 30mm thread, Blunt Point
50-1076-113B	Bone Screw, 6.0mm/110mm, 30mm thread, Blunt Point
50-1076-114B	Bone Screw, 6.0mm/110mm, 40mm thread, Blunt Point
50-1076-133B	Bone Screw, 6.0mm/130mm, 30mm thread, Blunt Point
50-1076-144B	Bone Screw, 6.0mm/140mm, 40mm thread, Blunt Point
50-1076-153B	Bone Screw, 6.0mm/150mm, 30mm thread, Blunt

